

AVAILABLE MATERIALS

SOFT MATERIALS

Material	Water content	UV Protection	Handling tint	Non-Ionic
SILICONE HYDROGEL	74%	NO	CLEAR	YES
SILICONE HYDROGEL	74%	YES	BLUE	YES
HYDROGEL	38%	AVAILABLE	CLEAR	YES
HYDROGEL	55%	NO	BLUE	NO
HYDROGEL	58%	AVAILABLE	CLEAR	YES
HYDROGEL	67%	AVAILABLE	CLEAR	YES
HYDROGEL	72%	NO	CLEAR	YES
HYDROGEL	77%	YES	CLEAR	YES
BENZ G3X	49%	YES	LIGHT BLUE	YES
CONTAMAC GM3	58%	YES	CLEAR	YES
CONTAMAC GM3	70%	YES	DARK BLUE	YES

FOR FURTHER DETAILS
CALL **0800 585115**
(FREEPHONE) OR EMAIL
ORDERS@ULTRAVISION.CO.UK



GAS PERMEABLE MATERIALS

Material type	Tint	Dk	Description / applications
First Choice Range identified by bold text			
OPTIMUM COMFORT	GREY	65	DAILY WEAR WITH ENHANCED WETTABILITY. WETTING ANGLE OF 6° FOR PROMOTION OF HIGH COMFORT LEVELS
OPTIMUM COMFORT	BLUE	65	
OPTIMUM COMFORT	GREEN	65	
OPTIMUM EXTRA	GREEN	100	AS FOR OPTIMUM COMFORT, BUT WITH HIGHER DK AND LOWER WETTING ANGLE 3°
OPTIMUM EXTRA	DARK BLUE	100	
OPTIMUM EXTREME	BLUE	125	AS FOR OPTIMUM COMFORT, BUT WITH HIGHER DK. WETTING ANGLE OF 6°
BOSTON EO	BLUE	58	MID-DK MATERIAL SUITABLE FOR DAILY WEAR APPLICATIONS.
BOSTON EO	BLUE UV	58	
BOSTON XO	VIOLET UV	100	
BOSTON IV	BLUE UV	19	WELL-PROVEN MATERIAL FOR DAILY WEAR APPLICATIONS.
BOSTON IV	CLEAR	19	
EQUALENS	BLUE	47	MID/HIGH DK MATERIAL FOR DAILY WEAR.
PF 45	BLUE	50	DEPENDABLE MID-DK GENERAL-PURPOSE MATERIAL FOR DAILY WEAR APPLICATIONS
PF 45	GREY	50	
HYBRID FS	BLUE	23	RGP CORE MATERIAL WITH A UNIQUE HYDROPHILIC SURFACE PROVIDING A WETTING ANGLE OF VIRTUALLY ZERO. SURFACES ARE NOT TREATED TO CREATE THIS EFFECT
HYBRID FS	GREY	23	
HYBRID FS PLUS	BLUE	60	SAME AS HYBRID FS BUT WITH HIGHER DK
HYBRID FS PLUS	GREY	60	
PARAGON® HDS 60	BLUE	60	MID-DK MATERIAL
PARAGON® HDS 60	GREEN	60	
PARAGON® HDS 100	BLUE	100	HIGH-DK MATERIAL
FLUOROPERM® 30	GREY	30	
FLUOROPERM® 60	BLUE	43	
AL02	LIGHT BLUE/CLEAR	28	
AL03	MID BLUE/CLEAR	50	

Note: DK measurement units ($\times 10^{-11}$ (cm²/sec)[ml O₂/(ml x mm Hg)])



UltraVision CLPL

Commerce Way, Leighton Buzzard, Bedfordshire LU74RW, United Kingdom

UK Orderline 0800 585115 (Freephone)

t +44 (0) 1525 381112 f +44 (0) 1525 370091
info@ultravision.co.uk ultravision.co.uk



© 2016 UltraVision International Limited
Registered in England and Wales no. 1408851 VAT No: GB 2519252 10
Definitive™ is a trademark of Contamac Limited