

# COSMETIC & PROSTHETIC

GIVING PEOPLE THE CONFIDENCE TO FACE THE WORLD

A WIDE RANGE OF TINT OPTIONS AND ENHANCEMENTS  
TO RESTORE THE NATURAL LOOK OF THE EYE



OUR INNOVATIVE DESIGNS MANUFACTURED FROM  
AN EXTENSIVE RANGE OF MATERIALS

## FEATURES

THREE MATERIAL OPTIONS

SIMULATED IRIDES

COLLARETTES AND  
LIMBAL BANDS

TINT DENSITIES:  
5% 10% 20% 30% 40%

## AVAILABLE LENS DESIGNS

DURAWAVE


IGEL

PRIMA


PRESTO

KERASOFT IC 77%


# COSMETIC & PROSTHETIC




SIMULATED IRIS S1




SIMULATED IRIS S1 WITH LIMBAL BAND



SIMULATED IRIS S2



SIMULATED IRIS S2 WITH COLLARETTE



SIMULATED IRIS S2 WITH COLLARETTE + LIMBAL BAND

## OPTIONAL WARRANTY

90 day warranty includes 2 exchanges and 75% Credit.

## EASY-FIT SYSTEM

Just call Customer Services with:

SPECTACLE Rx

K READINGS

HVID

BVD

and UltraVision will calculate the lens for you.

## TRIAL SETS

Clinical Services will be glad to talk you through the options available.

## COLOURS

Brown, Green Brown, Grey, Blue Grey, Blue, Aqua, Green, Yellow, Orange, Red, Violet, Black.

## PUPILS

Clear and Black pupils from 2.00 - 9.50 mm  
Opaque backing available.

PLEASE REQUEST A FULL PRODUCT BROCHURE FOR **LENS COLOUR CHART AND ORDERING INFORMATION**

FOR FURTHER DETAILS  
CALL **0800 585115**  
(FREEPHONE) OR EMAIL  
**ORDERS@ULTRAVISION.CO.UK**



**ultravision** CLPL  
INNOVATION IN PRACTICE

### UltraVision CLPL

Commerce Way, Leighton Buzzard, Bedfordshire LU74RW, United Kingdom

UK Orderline **0800 585115** (Freephone)

t +44 (0)1525 381112 f +44 (0)1525 370091  
info@ultravision.co.uk ultravision.co.uk



©2015 UltraVision International Limited  
Registered in England and Wales no. 1408851 VAT No: GB 322 5509 81