

TOPFLEX & EUREKA

CLASSIC RGP DESIGNS FROM ULTRAVISION



OUR INNOVATIVE DESIGNS MANUFACTURED FROM
AN EXTENSIVE RANGE OF MATERIALS

EUREKA

ENSURES EXCELLENT
FIRST TIME FITS

AVAILABLE AS SPHERES AND BACK
SURFACE TORIC DESIGN

TOPFLEX

DESIGN MINIMISES OVERALL
THICKNESS

AVAILABLE AS SPHERES AND FRONT
SURFACE TORIC DESIGN

FIT OPTIONS

12 LENS TRIAL SETS

EMPIRICAL ORDER SYSTEM

TOPFLEX & EUREKA

MATERIAL OPTIONS

OPTIMUM COMFORT	65DK
OPTIMUM EXTRA	100 DK
OPTIMUM EXTREME	125DK
SEE RGP MATERIAL LIST FOR FURTHER OPTIONS	

TOPFLEX & EUREKA PARAMETERS

BASE CURVES	7.00 TO 9.00 mm (0.05 steps)
DIAMETERS	8.50 TO 10.10 mm (0.30 steps)
SPHERICAL POWERS	+/- 30.00 DS (0.25 steps)
CYLINDER	-0.50 DC TO -6.00 DC
AXES	1° TO 180°

TRIAL SETS

12 lens Sphere and Toric* trial sets available with diameters of 9.50mm and base curves of 7.30 to 8.40 (in 0.10 steps).

*Back surface toric – Topflex. Front surface toric – Eureka

EMPIRICAL FITTING

Just call Customer Services with:

K READINGS

SPECTACLE REFRACTION

HVID + BVD

and UltraVision will calculate the lens for you.

OPTIONAL WARRANTY

Achieve the ideal fit without compromise and order a first fit lens with UltraVision's Fitting Warranty. Includes 2 free exchanges and a 75% credit for 90 days from the date of the original order.

ONLINE ORDERING

You can easily order replacement lenses through our online ordering system at ultravisionorders.co.uk

FOR FURTHER DETAILS
CALL **0800 585115**
(FREEPHONE) OR EMAIL
ORDERS@ULTRAVISION.CO.UK



ultravision CLPL
INNOVATION IN PRACTICE

UltraVision CLPL
Commerce Way, Leighton Buzzard, Bedfordshire LU74RW, United Kingdom

UK Orderline 0800 585115 (Freephone)

t +44 (0)1525 381112 f +44 (0)1525 370091
info@ultravision.co.uk ultravision.co.uk



©2015 UltraVision International Limited
Registered in England and Wales no. 1408851 VAT No: GB 322 5509 81