

# PRIMA 67

PROVIDE ALL YOUR PATIENT PESCRIPTION NEEDS  
WITH AN EASY-FIT, COST EFFECTIVE LENS



A MID-WATER CONTENT, COST-EFFECTIVE AND VERSATILE LENS

## PARAMETER RANGE

BASE CURVE 8.30/8.70/9.10 mm

DIAMETER 14.50 mm

SPHERES TO +/- 30.00 DS

## DESIGN FEATURES

HIGH WATER CONTENT, SPHERICAL TRI-CURVE  
FRONT SURFACE, SINGLE VISION CONTACT LENSES

# PRIMA 67

## MATERIAL

FILCON II	67% water, Filcon II 2
HANDLING TINT	Clear
DK	30 x 10-11 (cm <sup>2</sup> /sec) [ml O <sub>2</sub> /(ml x mm Hg)]

## MODALITY

Annual lens for daily wear only.

Pack Size Single Lenses

## OPTIONAL WARRANTY

Achieve a guaranteed fit with UltraVision's 90 day Fitting Warranty. This includes 2 free exchanges and a 75% credit from the date of the original order.

## ENCORE REPLACEMENT SCHEME

Once the contact lens fit is achieving the desired vision and comfort levels, patient details can be added to UltraVision's Encore™ Replacement Scheme which ensures a regular delivery of lenses\*, either to the practice or direct to patient. Registration with Encore™ ensures wear modality compliance and helps to encourage practitioner loyalty.

\*Please note: A letter confirming the prescription is sent to the practice at least six weeks prior to the manufacturing date of the lens. Lens prescriptions should be checked prior to the manufacturing date in order that the correct lenses are despatched. Alterations made to the prescription after the manufacturing date carries a £10 administration charge.

## ONLINE ORDERING

You can easily order replacement lenses through our online ordering system at [ultravisionorders.co.uk](http://ultravisionorders.co.uk)

FOR FURTHER DETAILS  
CALL **0800 585115**  
(FREEPHONE) OR EMAIL  
**ORDERS@ULTRAVISION.CO.UK**



**ultravision** CLPL  
INNOVATION IN PRACTICE

### UltraVision CLPL

Commerce Way, Leighton Buzzard, Bedfordshire LU74RW, United Kingdom

UK Orderline 0800 585115 (Freephone)

[info@ultravision.co.uk](mailto:info@ultravision.co.uk) [www.ultravision.co.uk](http://www.ultravision.co.uk)

©2015 Ultra Vision International Limited  
Registered in England and Wales no. 1408851 VAT No: GB 322 5509 81