

FULL PRODUCT LISTING

SOFT CONTACT LENSES FOR KERATOCONUS AND IRREGULAR CORNEAS

Design ¹	BCOR	Diameter ²	Power	Cyl	Add	Material	Modality
KERASOFT THIN	7.40 TO 9.40 (0.20 STEPS)	14.00 TO 15.50	+/- 30.00	-0.50 TO -15.00		SiH 74% & HYDROGEL 77%	3/12 MONTHS
KERASOFT IC	7.40 TO 9.40 (0.20 STEPS)	14.00 TO 15.50	+/- 30.00	-0.50 TO -15.00		SiH 74% & HYDROGEL 77%	3/12 MONTHS
KERASOFT 3	8.00 TO 8.60 (0.20 STEPS)	14.00/14.50/15.00	+/- 30.00	-0.50 TO -11.00	+1.00 TO +3.50	SiH 74%	3 MONTHS
KERASOFT 2	8.00 TO 8.60 (0.20 STEPS)	14.00/14.50/15.00	+/- 30.00	-0.50 TO -11.00	+1.00 TO +3.50	HYDROGEL 49%	UP TO 12 MONTHS

SOFT CONTACT LENSES (SPHERES, TORICS AND MULTIFOCALS)

HYDROWAVE SiH	7.00 TO 9.60 (0.10 STEPS)	12.50 TO 16.00	+/- 30.00	-0.50 TO -11.00	+1.00 TO +3.50	SiH 74%	3 MONTHS
DURAWAVE	7.80 TO 9.40 (0.10 STEPS)	13.00 TO 16.00	+/- 30.00	-0.50 TO -11.00	+1.00 TO +3.50	HYDROGEL 38%, 49%, 55%, 58%, 77% & GM3	UP TO 12 MONTHS
SAM MF SOFT	7.00 TO 9.40 (0.10 STEPS)	12.50 TO 16.00	+/- 30.00	-0.50 TO -6.00	+1.00 TO +3.50	HYDROGEL 58% & 72%	UP TO 12 MONTHS
IGEL	7.00 TO 9.40 (0.10 STEPS)	12.50 TO 16.00	+/- 30.00	-0.50 TO -11.00		HYDROGEL 38%, 49%, 58%, 67% & 77%	UP TO 12 MONTHS
SIMPLON	6.70 TO 9.40 (0.10 STEPS)	11.00 TO 16.50	+/- 30.00	-0.50 TO -11.00		HYDROGEL 38%, 72%, 77% & GM3	UP TO 12 MONTHS
RXT	7.00 TO 9.40 (0.10 STEPS)	12.50 TO 14.50	+20.00 TO -30.00	-0.50 TO -11.00		HYDROGEL 38%	UP TO 12 MONTHS
STANDARD BANDAGE	8.30/8.70/9.10	14.00/14.50/15.00	AFOCAL			SiH 74% HYDROGELS: 67%, 72% & 77%	3/12 MONTHS
CUSTOM BANDAGE	7.00 TO 9.40 (0.10 STEPS)	12.50 TO 19.00	AFOCAL			SiH 74% HYDROGELS: 67%, 72% & 77%	3/12 MONTHS
PAEDIATRIC SiH	6.50 TO 8.20 (0.10 STEPS)	11.50 TO 14.00	+/- 45.00			SiH 74%	3 MONTHS
PAEDIATRIC	6.50 TO 8.20 (0.10 STEPS)	11.50 TO 14.00	+/- 45.00			HYDROGEL 67%, 72%, 77% & GM3	UP TO 12 MONTHS
APHAKIC & MYOPIC SiH	7.00 TO 9.80 (0.10 STEPS)	12.50 TO 16.00	+/- 45.00			SiH 74%	3 MONTHS
CD 38	8.30/8.50/8.70/8.90	14.00	+/- 20.00			HYDROGEL 38%	UP TO 12 MONTHS
PRESTO 58	8.40/8.80	14.50	+/- 20.00			HYDROGEL 58%	UP TO 12 MONTHS
PRIMA 67	8.30/8.70/9.10	14.50	+/- 20.00			HYDROGEL 67%	UP TO 12 MONTHS
MINI (SOFTINT ALSO AVAILABLE)	7.90 TO 9.00 (0.10 STEPS)	12.50/13.00/13.50	+/- 30.00			HYDROGEL 38%	UP TO 12 MONTH
H67 (SOFTINT ALSO AVAILABLE)	8.00 TO 9.40 (0.10 STEPS)	14.00	+/- 30.00			HYDROGEL 67%	UP TO 12 MONTHS
OMNIFLEX	8.10 TO 9.40 (0.10 STEPS)	14.30/14.50	+/- 20.00	-0.50 TO -6.00		HYDROGEL 67%	UP TO 12 MONTHS
ZERO 6 SOFTINTS	8.10 TO 9.40 (0.10 STEPS)	14.00/14.50	+/- 25.00	-0.50 TO -6.00		HYDROGEL 38%	UP TO 12 MONTHS
OMEGA 38	8.60	14.00	+6.00 TO -12.00			HYDROGEL 38%	UP TO 6 MONTHS
U55 AB	8.70	14.50	+6.00 TO -12.00			HYDROGEL 55%	UP TO 3 MONTHS

FULL PRODUCT LISTING

RGP CONTACT LENSES ³

Design ¹	BCOR	Diameter ²	Power	Cyl	Add	Material	Modality
EUREKA	7.00 TO 9.00 (0.05 STEPS)	8.50 TO 10.10	+/- 30.00	-0.50 TO -6.00		OPTIMUM COMFORT/EXTRA/EXTREME	UP TO 12 MONTHS
SAM MF RGP	7.00 TO 9.00 (0.05 STEPS)	8.90 TO 10.10	+/- 30.00	-0.50 TO -6.00	+1.00 TO +4.50	OPTIMUM COMFORT/EXTRA/EXTREME	UP TO 12 MONTHS
TOPFLEX	7.00 TO 9.00 (0.05 STEPS)	8.50 TO 10.10	+/- 30.00	-0.50 TO -6.00		OPTIMUM COMFORT/EXTRA/EXTREME	UP TO 12 MONTHS
XTRALENS	7.00 TO 9.00 (0.05 STEPS)	9.00 TO 11.70 (0.30 STEPS)	+/- 30.00			OPTIMUM COMFORT/EXTRA/EXTREME	UP TO 12 MONTHS

RGP CONTACT LENSES FOR KERATOCONUS ³

ALK	5.50 TO 8.20 (0.05 STEPS)	8.50 TO 11.00	+/- 30.00	-0.50 TO -6.00		OPTIMUM COMFORT/EXTRA/EXTREME	UP TO 12 MONTHS
WOODWARD	5.00 TO 8.60 (0.05 STEPS)	8.50 TO 11.00	+/- 30.00	-0.50 TO -6.00		OPTIMUM COMFORT/EXTRA/EXTREME	UP TO 12 MONTHS
MCGUIRE	5.00 TO 8.60 (0.05 STEPS)	8.50 TO 11.00	+/- 30.00	-0.50 TO -6.00		OPTIMUM COMFORT/EXTRA/EXTREME	UP TO 12 MONTHS

MONTHLY DISPOSABLE

ULTRAWAVE	8.70	14.40	+6.00 TO -12.00			HYDROGEL 57%	1 MONTH
ULTRAWAVE T	8.70	14.50	+4.00 TO -8.00	-1.00 / -1.75		HYDROGEL 57%	1 MONTH
SPECIALTY 55 AB	8.70	14.50	+6.00 TO -12.00			HYDROGEL 55%	1 MONTH

1. Available parameters are material dependant

2. All lens diameters are available in 0.10 steps unless specified

3. Other RGP materials are available

If you require any parameters outside of the listed range, please contact us

UltraVision CLPL

Commerce Way, Leighton Buzzard, Bedfordshire LU7 4RW, United Kingdom

UK Orderline 0800 585115 (Freephone)

t +44 (0) 1525 381112 f +44 (0) 1525 370091
info@ultravision.co.uk ultravision.co.uk



©2016 UltraVision International Limited
Registered in England and Wales no. 1408851 VAT No: GB 322 5509 81