



# MAKE SENSE OF KERATOCONUS

KERA  SOFT THIN

ULTRA-THIN SPECIALIST CONTACT LENS FOR KERATOCONUS AND IRREGULAR CORNEAS

# KERASOFT® THIN ALLOW YOUR PATIENTS TO ENJOY LIFE TO THE FULL

For patients with keratoconus, lifestyle often takes second place to ensuring they have the visual acuity necessary for work or driving. Rigid contact lens options can mean a compromise on comfort, leading to reduced wear time which generally has an impact on leisure activities.

KeraSoft® Thin delivers visual acuity comparable to that of rigid options with all day comfort that ensures all day wear, every day, whatever your patient wants to do.

**BE KIND TO YOUR PATIENTS' EYES WITH KERASOFT® THIN**  
ULTRA-THIN SPECIALIST CONTACT LENS FOR KERATOCONUS AND IRREGULAR CORNEAS

## FOR ALL TYPES OF KERATOCONUS & REVERSE GEOMETRY CORNEAS

### ULTRA-THIN DESIGN DELIVERS ENHANCED PERIPHERAL FIT

KeraSoft® Thin's unique, thinner design increases lens flexibility, allowing it to drape better over the peripheral cornea. This directly reduces the need for peripheral changes on complex corneas such as Post Graft and Nipple cones. Compared to thicker soft lens designs, KeraSoft® Thin is quicker and easier to fit, thus reducing chair time.

### COMPLEX CASES BENEFIT FROM UNIQUE KERASOFT® THIN DESIGN

All your patients can benefit from the unique design of the KeraSoft® Thin lens. Over twenty years of research and development has proven that fitting thin can provide optimal VA together with lifestyle enhancing comfort for even very complex corneas. Improved drapage reduces tear pooling under the lens which in turn stabilises visual acuity.

### THINNER LENSES DELIVER INCREASED OXYGEN TO THE CORNEA

Increased oxygen supply helps maintain a healthy cornea - particularly important for post graft cases. Many cases of keratoconus require high minus powers which adds overall thickness to the final powered lens. KeraSoft® Thin works to reduce this thickness and therefore improves wearing time.



# BE KIND TO KERATOCONUS

## FIT SOFT | FIT HEALTHY | FIT THIN

**KERASOFT<sup>®</sup> THIN** ULTRA-THIN SPECIALIST CONTACT LENS FOR KERATOCONUS AND IRREGULAR CORNEAS

**KERASOFT<sup>®</sup> THIN**  
**EASY TO FIT**  
**OPTIMAL COMFORT**  
**HEALTHIER EYES**

Eight Lens Fitting Set  
First Choice Lens Calculator

The MoRoCCo VA fitting method guides you to achieve optimal fit in minimum chair time using the eight Lens Fitting Set. Additionally, First Choice Lens can be calculated from a Topography Map - please refer to the Fitting Guide for more details.

### PARAMETER RANGE

Base Curve	7.40 TO 9.40 mm
Diameter	14.00 14.50 15.00mm
Sphere	+/- 30.00 DS
Cyl	-0.50 to -15.00 DC
Axis	1°- 180° (1° steps)
Peripheries	FLT4 to FLT1 STD STP1 to STP4
SMC*	Up to 2 sectors

### MATERIAL


Filcon V 3 SiH, 74% water\*\* (3 month replacement)  
Filcon II 3, 77% Water (12 month replacement)

\*Sector Management Control \*\*Definitive™ SiH

# ULTRA-THIN DESIGN INCREASES OXYGEN PERMEABILITY & COMFORT

FOR FURTHER DETAILS CALL OUR FREEPHONE UK ORDERLINE

---

 **0800 585115**

---

Tel +44 (0) 1525 381112 • [info@ultravision.co.uk](mailto:info@ultravision.co.uk) • [www.ultravision.co.uk](http://www.ultravision.co.uk)

**KERA**  **SOFT** **THIN**

ULTRA-THIN SPECIALIST CONTACT LENS FOR KERATOCONUS AND IRREGULAR CORNEAS