

# McGUIRE

## KERATOCONIC RGP DESIGN FROM ULTRAVISION



OUR INNOVATIVE DESIGNS MANUFACTURED FROM AN EXTENSIVE RANGE OF MATERIALS

THREE POINT TOUCH FITTING FOR KERATOCONUS

DESIGN BASED ON THE RELATIONSHIP BETWEEN THE RADII OF CONE AND RATE OF FLATTENING AT THE ADJOINING AND PERIPHERAL CORNEAL AREAS

# McGUIRE

## MATERIAL OPTIONS\*

OPTIMUM COMFORT	65 DK
OPTIMUM EXTRA	100 DK
OPTIMUM EXTREME	125 DK

\*Other RGP materials available on request

## PARAMETER RANGE\*

BASE CURVES	5.00 TO 8.60 mm (0.05 STEPS)
DIAMETERS	8.50 TO 11.00 mm (0.05 STEPS)
POWERS	- 30.00 TO + 30.00 DC

\*Toric and multifocal designs available on request

## EASY TO ORDER FOR EMPIRICAL FITS

Just call Customer Services with:

SPECTACLE Rx	K READINGS
HVID	BVD

and UltraVision will calculate the lens for you

## TRIAL SETS

Early, moderate and advanced fitting sets available

## MODALITY

Single lens daily wear, up to 12 months.

## ULTRAVISION'S ENCORE™ REPLACEMENT SCHEME

Once the contact lens fit is achieving the desired vision and comfort levels, patient details can be added to UltraVision's Encore™ Replacement Scheme which ensures a regular delivery of lenses\*, either to the practice or direct to patient. Registration with Encore™ ensures wear modality compliance and helps to encourage practitioner loyalty.

\*Please note: A letter confirming the prescription is sent to the practice at least six weeks prior to the manufacturing date of the lens. Lens prescriptions should be checked prior to the manufacturing date in order that the correct lenses are despatched. Alterations made to the prescription after the manufacturing date carries a £10 administration charge.

## OPTIONAL WARRANTY

Achieve a guaranteed fit with UltraVision's 90 day Fitting Warranty. This includes 2 free exchanges and a 75% credit from the date of the original order.

## FITTING MCGUIRE FROM A TRIAL SET

Perform over refraction on chosen fitting lens and call Customer Services with details. UltraVision will calculate the lens for you. Fitting guides are available at [ultravision.co.uk](http://ultravision.co.uk).

FOR FURTHER DETAILS  
CALL **0800 585115**  
(FREEPHONE) OR EMAIL  
**ORDERS@ULTRAVISION.CO.UK**



UltraVision CLPL  
Commerce Way, Leighton Buzzard, Bedfordshire LU74RW, United Kingdom

UK Orderline 0800 585115 (Freephone)  
[info@ultravision.co.uk](mailto:info@ultravision.co.uk) [ultravision.co.uk](http://ultravision.co.uk)



©2015 UltraVision International Limited  
Registered in England and Wales no. 1408851 VAT No: GB 322 5509 81