

# PAEDIATRIC

DESIGNED IN CONJUNCTION WITH A LEADING NHS HOSPITAL

HIGH WATER CONTENT CONTACT LENSES IN DESIGNS SPECIFICALLY FOR THE VERY YOUNG EYE



OUR INNOVATIVE DESIGNS MANUFACTURED FROM AN EXTENSIVE RANGE OF MATERIALS

## TWO MATERIAL OPTIONS

LATHE-CUT SiH, 74%  
WATER CONTENT

FILCON II 2,  
77% WATER CONTENT\*

\*Other materials available

## APPLICATIONS

PAEDIATRIC APHAKIA

HIGH ANISOMETROPIA

HIGH REFRACTIVE ERRORS

INABILITY TO  
WEAR SPECTACLES

## FITTING GUIDE

Written by Jon Whittle BSc(Hons) MSc PhD,  
Lead Optometrist Sheffield  
Children's Teaching Hospital / Royal  
Hallamshire Hospital, Sheffield.

Available at [ultravision.co.uk](http://ultravision.co.uk)

# PAEDIATRIC

## MODALITY

SiH - Daily wear, 3 month lens.  
Filcon II 2 - Daily wear, 12 month lens

## AVAILABILITY

Available as a single lens.  
SiH also supplied as a cost effective 2-pack option

## READINGS UNAVAILABLE?

If you cannot manage to obtain any readings, the following can be used as a guide for a first fit lens:

AGE	BCOR	DIAMETER	POWER
0-2 MONTHS	7.00 mm	11.50 mm	+43.00 DS
2-4 MONTHS	7.30 mm	12.00 mm	+38.00 DS
4-8 MONTHS	7.50 mm	12.50 mm	+34.00 DS
8-12 MONTHS	7.70 mm	13.00 mm	+30.00 DS

## PARAMETER RANGE

BASE CURVE	6.50 TO 8.20 mm
DIAMETER	11.50 TO 14.00 mm
POWERS	+/- 45.00 DS

## EASY-FIT SYSTEM

Just call Customer Services with:

K READINGS

HVID

REFRACTION

and UltraVision will calculate the lens for you.

## WARRANTY

SiH Paediatric automatically includes fitting warranty.  
For non SiH materials, an optional warranty can be purchased with the lens.

2-packs carry manufacturing warranty only

## ONLINE ORDERING

You can easily order replacement lenses through our online ordering system at [ultravisionorders.co.uk](http://ultravisionorders.co.uk)

FOR FURTHER DETAILS  
CALL **0800 585115**  
(FREEPHONE) OR EMAIL  
**ORDERS@ULTRAVISION.CO.UK**



**UltraVision CLPL**  
Commerce Way, Leighton Buzzard, Bedfordshire LU74RW, United Kingdom

**UK Orderline 0800 585115** (Freephone)  
[info@ultravision.co.uk](mailto:info@ultravision.co.uk) [ultravision.co.uk](http://ultravision.co.uk)



©2015 UltraVision International Limited  
Registered in England and Wales no. 1408851 VAT No: GB 322 5509 81