

IGEL™

DAILY WEAR CONTACT LENSES

PRECISION LATHE CONTACT LENSES



VISION CORRECTION WITH MAXIMUM STABILITY AND COMFORT

PARAMETER RANGE

BASE CURVE	7.00mm to 9.40mm
DIAMETER	12.50mm to 16.00mm
POWER RANGE	+/-30.00DS
CYLINDER RANGE	-0.50 to -11.00
AXIS	1 to 180

LENS DESIGN

Igel is available in both toric and spherical lenses; designed for patients who cannot tolerate aberration control.

Optical prism is available with the toric lens design. Advantages of adding optical prism include; correcting Diplopia (double vision) and some muscle imbalances, reducing headaches.

IGEL™

MATERIAL OPTIONS*

SELECT BENZ (HIOXIFILCON B1), 49% WATER CONTENT

GM3, 58% CLEAR, 72% BLUE WATER CONTENT

METHAFILCON 55% WATER CONTENT

HYDROGEL (FILCON II 2), 38%, 58%, 67%, 72%, 77% WATER CONTENT

*Igel™ design can be added to any soft lens material on request.

EASY-FIT SYSTEM - NO TRIAL SET NEEDED

Just call Customer Services with:

SPECTACLE Rx

K READINGS

HVID

BVD

and UltraVision will calculate the lens for you.

MODALITY

Single lens daily wear, up to 12 months.

ENCORE REPLACEMENT SCHEME

Ensure a regular delivery of lenses, either to the practice or direct to patient. Registration with Encore ensures wear modality compliance and helps retain practitioner loyalty.

Simply ask our Customer Service team to add the patient details to UltraVision's Encore replacement scheme.

ONLINE ORDERING

You can easily order replacement lenses through our online ordering system at ultravisionorders.co.uk

OPTIONAL WARRANTY

Achieve the ideal fit without compromise and order a first fit lens with UltraVision's Fitting Warranty. Includes 2 free exchanges and a 75% credit for 90 days from the date of the original order.

FOR FURTHER DETAILS
CALL **0800 585115 (FREEPHONE)**

OR EMAIL
ORDERS@ULTRAVISION.CO.UK



UltraVision CLPL

Commerce Way, Leighton Buzzard, Bedfordshire, LU7 4RW, United Kingdom

UK Orderline 0800 585115 (Freephone)

t+44 (0)1525 381112 f+44 (0)1525 370091
info@ultravision.co.uk www.ultravision.co.uk



©2015 UltraVision International Limited
Registered in England and Wales no. 1408851 VAT No: GB 322 5509 81

M1080/0514/REV. 2

