

BANDAGE™

Protection & Comfort for Damaged or Compromised Corneas

ASSISTS HEALING

- Use after corneal trauma including, but not limited to;
- Post-Surgery
 - Stabilisation of Epithelium
 - Bullous Keratopathy
 - Mechanical Protection
 - Stevens-Johnson Syndrome
 - Aqueous Leaks
 - Keratoplasty

VISUAL ACUITY IMPROVEMENTS

- Thygeson's Superficial Punctuate Keratitis

STANDARD RANGE*

Base Curve 8.30 / 8.70 / 9.10 mm
Diameters 14.00 / 14.50 / 15.00 mm

PAIN RELIEF

- Bullous Keratopathy
- Fuchs' Dystrophy
- Indolent Corneal Ulcer*
- Exposure via Trauma/Injury
- Entropion
- Trichiasis

SPLINT

- Corneal Injury
- Corneal Thinning
- Corneal Perforation

CUSTOM RANGE

Base Curve 7.00 to 9.40 mm
Diameters 12.50 to 19.00 mm

PHYSICAL BARRIER

- Filamentary Keratitis
- Chemical/Thermal burns
- Trachoma
- Stevens-Johnson syndrome

POST-OCULAR SURGERY

- Leaking blebs
- Pain relief after corneal surgical techniques
- Protection from suture ends
- Recurrent corneal erosion syndrome
- Anterior Membrane Dystrophies

MATERIAL OPTIONS

Definitive® SiH, 74% water content
Filcon II 2, 77% water content
Filcon II 2, 67% water content

*Standard range Bandage™ lenses are only available in Filcon II 2, 77% water content and Filcon II 2, 67% water content
Bandage lenses are designed for single use only

FOR FURTHER DETAILS
FREEPHONE **0800 585 115**
CALL **+44 (0) 1525 381112**
EMAIL **ORDERS@ULTRAVISION.CO.UK**
WEBSITE **ULTRAVISION.CO.UK**

