

ULTRAWAVE[®]

SUPERB COMFORT AND ENHANCED VISION

BIOMIMETIC MATERIAL AND UNIQUE ASPHERIC DESIGN IN A MONTHLY DISPOSABLE LENS



OUR INNOVATIVE DESIGNS MANUFACTURED FROM AN EXTENSIVE RANGE OF MATERIALS

FEATURES

BIOMIMETIC MATERIAL
ENSURES SUPERB
ALL DAY COMFORT

VISUAL BENEFITS
ENHANCED CONTRAST
SENSITIVITY INCREASED
DEPTH PERCEPTION GLARE
REDUCTION

PRESBYOPIA CORRECTION

ASPHERIC DESIGN WITH UP
TO +1.00DS NEAR ADD EFFECT

ULTRAWAVE[®]
SPHERICAL POWERS MASK
UP TO 1.25 DC


ULTRAWAVE[®] TORIC
ASPHERIC DESIGN ENHANCES
ASTIGMATIC CORRECTION

EASY-FIT

DIAGNOSTIC LENSES
AND FITTING
BANKS AVAILABLE

ULTRAWAVE®

SPHERICAL ABERRATION CORRECTION



UltraVision's patented aberration control system incorporates the results of studies of over one thousand eyes to arrive at a design that minimises the blurring effect from spherical aberration.

MATERIAL

Biomimetic for greater comfort

43% HIOXIFILCON 57% WATER CONTENT

MODALITY

Daily wear, 1 month lenses

PARAMETER RANGE

ULTRAWAVE®

BASE CURVE/DIAMETER

8.70 : 14.40 mm

SPHERICAL POWERS

-12.00 TO +6.00DS

ULTRAWAVE® TORIC

BASE CURVE/DIAMETER

8.70 : 14.50 mm

SPHERICAL POWERS

-8.00 TO +4.00DS

CYLINDER POWERS

-1.00 & -1.75 DC

AXES

10, 20, 80, 90, 100, 160, 170, 180

ENCORE REPLACEMENT SCHEME

Ensure a regular delivery of lenses, either to the practice or direct to patient. Registration with Encore ensures wear modality compliance and helps retain practitioner loyalty.

Simply ask our Customer Service team to add the patient details to UltraVision's Encore replacement scheme.

ONLINE ORDERING

You can easily order replacement lenses through our online ordering system at ultravisionorders.co.uk

CONTACT US FOR **FREE**
DIAGNOSTIC LENSES

FOR FURTHER DETAILS
CALL **0800 585115**
(FREEPHONE) OR EMAIL
ORDERS@ULTRAVISION.CO.UK



ultravision CLPL
INNOVATION IN PRACTICE

UltraVision CLPL

Commerce Way, Leighton Buzzard, Bedfordshire LU74RW, United Kingdom

UK Orderline 0800 585115 (Freephone)

info@ultravision.co.uk ultravision.co.uk

©2015 UltraVision International Limited
Registered in England and Wales no. 1408851 VAT No: GB 322 5509 81

M1038/1112/REV.2