

FULL PRODUCT LISTING

SOFT CONTACT LENSES FOR KERATOCONUS AND IRREGULAR CORNEAS

Design ¹	BCOR	Diameter ²	Power	Cyl	Add	Material	Modality
KERASOFT THIN	7.40 TO 9.40 (0.20 STEPS)	14.00 TO 15.50	+/- 30.00	-0.50 TO -15.00		SiH 74% CLEAR & HYDROGEL 77%	3/12 MONTHS
KERASOFT IC	7.40 TO 9.40 (0.20 STEPS)	14.00 TO 15.50	+/- 30.00	-0.50 TO -15.00		SiH 74% CLEAR & HYDROGEL 77%	3/12 MONTHS
KERASOFT 3	8.00 TO 8.60 (0.20 STEPS)	14.00/14.50/15.00	+/- 30.00	-0.50 TO -11.00	+1.00 TO +3.50	SiH 74% CLEAR	3 MONTHS
KERASOFT 2	8.00 TO 8.60 (0.20 STEPS)	14.00/14.50/15.00	+/- 30.00	-0.50 TO -11.00	+1.00 TO +3.50	HYDROGEL 49%	UP TO 12 MONTHS

SOFT CONTACT LENSES (SPHERES, TORICS AND MULTIFOCALS)

HYDROWAVE SiH	7.00 TO 9.60 (0.10 STEPS)	12.50 TO 16.00	+/- 30.00	-0.50 TO -11.00	+1.00 TO +3.50	SiH 74% Blue	3 MONTHS
DURAWAVE*	7.80 TO 9.40 (0.10 STEPS)	13.00 TO 16.00	+/- 30.00	-0.50 TO -11.00	+1.00 TO +3.50	HYDROGEL 38%, 49%, 55%, 58%, 77% & GM3 58	UP TO 12 MONTHS
SAM MF SOFT	7.00 TO 9.40 (0.10 STEPS)	12.50 TO 16.00	+/- 30.00	-0.50 TO -6.00	+1.00 TO +3.50	HYDROGEL 58% & 72%	UP TO 12 MONTHS
IGEL*	7.00 TO 9.40 (0.10 STEPS)	12.50 TO 16.00	+/- 30.00	-0.50 TO -11.00		HYDROGEL 38%, 49%, 58%, 67% & 77%	UP TO 12 MONTHS
SIMPLON	6.70 TO 9.40 (0.10 STEPS)	11.00 TO 16.50	+/- 30.00	-0.50 TO -11.00		HYDROGEL 38%, 72%, 77% & GM3 58	UP TO 12 MONTHS
RXT	7.00 TO 9.40 (0.10 STEPS)	12.50 TO 14.50	+20.00 TO -30.00	-0.50 TO -11.00		HYDROGEL 38%	UP TO 12 MONTHS
STANDARD BANDAGE	8.30/8.70/9.10	14.00/14.50/15.00	AFOCAL			SiH 74% HYDROGELS: 67%, 72% & 77%	3/12 MONTHS
CUSTOM BANDAGE	7.00 TO 9.40 (0.10 STEPS)	12.50 TO 19.00	AFOCAL			SiH 74% HYDROGELS: 67%, 72% & 77%	3/12 MONTHS
SPECIALIST BANDAGE	8.40/8.70/9.00	16.00/18.00	AFOCAL			HYDROGEL 77%	1 MONTH
PAEDIATRIC SiH	6.50 TO 8.20 (0.10 STEPS)	11.50 TO 14.00	+/- 45.00			SiH 74% Blue	3 MONTHS
PAEDIATRIC	6.50 TO 8.20 (0.10 STEPS)	11.50 TO 14.00	+/- 45.00			HYDROGEL 67%, 72%, 77% & GM3 58	UP TO 12 MONTHS
APHAIC & MYOPIC SiH	7.00 TO 9.80 (0.10 STEPS)	12.50 TO 16.00	+/- 45.00			SiH 74% Blue	3 MONTHS
CD 38	8.30/8.50/8.70/8.90	14.00	+/- 20.00			HYDROGEL 38%	UP TO 12 MONTHS
PRESTO 58 *	8.40/8.80	14.50	+/- 20.00			HYDROGEL 58%	UP TO 12 MONTHS
PRIMA 67 *	8.30/8.70/9.10	14.50	+/- 20.00			HYDROGEL 67%	UP TO 12 MONTHS
MINI (SOFTINT ALSO AVAILABLE)	7.90 TO 9.00 (0.10 STEPS)	12.50/13.00/13.50	+/- 30.00			HYDROGEL 38%	UP TO 12 MONTH
H67 (SOFTINT ALSO AVAILABLE)	8.00 TO 9.40 (0.10 STEPS)	14.00	+/- 30.00			HYDROGEL 67%	UP TO 12 MONTHS
OMNIFLEX SOFTINTS	8.10 TO 9.40 (0.10 STEPS)	14.30/14.50	+/- 20.00	-0.50 TO -6.00		HYDROGEL 67%	UP TO 12 MONTHS
ZERO 6 SOFTINTS	8.10 TO 9.40 (0.10 STEPS)	14.00/14.50	+/- 25.00	-0.50 TO -6.00		HYDROGEL 38%	UP TO 12 MONTHS

* Tint available upon request

FULL PRODUCT LISTING

RGP CONTACT LENSES

Design ¹	BCOR	Diameter ²	Power	Cyl	Add	Material	Modality
EUREKA	7.00 TO 9.00 (0.05 STEPS)	8.50 TO 10.10	+/- 30.00	-0.50 TO -6.00		AVAILABLE IN ALL RGP MATERIALS ³	UP TO 12 MONTHS
SAM MF RGP	7.00 TO 9.00 (0.05 STEPS)	8.90 TO 10.10	+/- 30.00	-0.50 TO -6.00	+1.00 TO +4.50	AVAILABLE IN ALL RGP MATERIALS ³	UP TO 12 MONTHS
TOPFLEX	7.00 TO 9.00 (0.05 STEPS)	8.50 TO 10.10	+/- 30.00	-0.50 TO -6.00		AVAILABLE IN ALL RGP MATERIALS ³	UP TO 12 MONTHS
XTRALENS	7.00 TO 9.00 (0.05 STEPS)	9.00 TO 11.70 (0.30 STEPS)	+/- 30.00			AVAILABLE IN ALL RGP MATERIALS ³	UP TO 12 MONTHS

RGP CONTACT LENSES FOR KERATOCONUS³

ALK	5.50 TO 8.20 (0.05 STEPS)	8.50 TO 11.00	+/- 30.00	-0.50 TO -6.00		AVAILABLE IN ALL RGP MATERIALS ³	UP TO 12 MONTHS
WOODWARD	5.00 TO 8.60 (0.05 STEPS)	8.50 TO 11.00	+/- 30.00	-0.50 TO -6.00		AVAILABLE IN ALL RGP MATERIALS ³	UP TO 12 MONTHS
MCGUIRE	5.00 TO 8.60 (0.05 STEPS)	8.50 TO 11.00	+/- 30.00	-0.50 TO -6.00		AVAILABLE IN ALL RGP MATERIALS ³	UP TO 12 MONTHS

MONTHLY DISPOSABLE

AVANTI	8.00 TO 9.20 (0.30 STEPS)	14.00 TO 15.00 (0.50 STEPS)	+/-20.00	-0.75 TO -8.00	+0.50 TO +4.00	SiH 74% BLUE	1 MONTH
ULTRAWAVE	8.60	14.20	+6.00 TO -12.00			HYDROGEL HIOXIFILCON 57%	1 MONTH
ULTRAWAVE T	8.60	14.20	+4.00 TO -8.00	-0.75 / -1.75		HYDROGEL HIOXIFILCON 57%	1 MONTH
SPECIALTY 55 AB	8.70	14.50	+6.00 TO -12.00			HYDROGEL 55%	1 MONTH

1. Available parameters are material dependant

2. All lens diameters are available in 0.10 steps unless specified

3. Please refer to UltraVision's Material List Document

If you require any parameters outside of the listed range, please contact us

UltraVision CLPL

Commerce Way, Leighton Buzzard, Bedfordshire LU7 4RW, United Kingdom

UK Orderline 0800 585115 (Freephone)

t +44 (0) 1525 381112 f +44 (0) 1525 370091
info@ultravision.co.uk ultravision.co.uk



©2018 UltraVision International Limited
Registered in England and Wales no.1408851 VAT No: GB 251 925210