

ALKTM

KERATOCONIC RGP DESIGN FROM ULTRAVISION



CUSTOM MADE LENSES WITH ABERRATION CONTROL FOR OPTIMAL VISION

LENS SPECIFIC EDGE DESIGN

ABERRATION FREE FRONT OPTICAL CURVES

CUSTOM MADE FOR THE KERATOCONIC EYE

BASED ON KERATOCONUS SEVERITY, NOT K-READING MEASUREMENTS

ALK™

MATERIAL OPTIONS*

OPTIMUM COMFORT	65 DK
OPTIMUM EXTRA	100 DK
OPTIMUM EXTREME	125 DK

*Other RGP materials available on request

PARAMETER RANGE

BASE CURVE	6.00 to 8.20mm
DIAMETER	8.50 to 11.00mm
POWER RANGE	+/-30.00DS

*Toric designs available on request

FITTING ALK

Fitting sets are available. An over refraction should be performed on the chosen fitting lens. The Customer Service team at UltraVision will calculate the lens for you. Fitting guides are available on ultravision.co.uk.

MODALITY

Single lens daily wear, up to 12 months.

ULTRAVISION'S ENCORE™ REPLACEMENT SCHEME

Once the contact lens fit is achieving the desired vision and comfort levels, patient details can be added to UltraVision's Encore™ Replacement Scheme which ensures a regular delivery of lenses*, either to the practice or direct to patient. Registration with Encore™ ensures wear modality compliance and helps to encourage practitioner loyalty.

*Please note: A letter confirming the prescription is sent to the practice at least six weeks prior to the manufacturing date of the lens. Lens prescriptions should be checked prior to the manufacturing date in order that the correct lenses are despatched. Alterations made to the prescription after the manufacturing date carries a £10 administration charge.

OPTIONAL WARRANTY

Achieve a guaranteed fit with UltraVision's 90 day Fitting Warranty. This includes 2 free exchanges and a 75% credit from the date of the original order.

FOR FURTHER DETAILS
CALL **0800 585115 (FREEPHONE)**

OR EMAIL
ORDERS@ULTRAVISION.CO.UK



UltraVision CLPL

Commerce Way, Leighton Buzzard, Bedfordshire, LU7 4RW, United Kingdom

UK Orderline 0800 585115 (Freephone)

t+44 (0)1525 381112 t+44 (0)1525 370091
info@ultravision.co.uk www.ultravision.co.uk



©2015 UltraVision International Limited
Registered in England and Wales no. 1408851 VAT No: GB 322 5509 81

MIO86/0814/REV.1

