

# COSMETIC & PROSTHETIC

## TINTED LENS COLOUR CHART

	BROWN	GREEN BROWN	GREY	BLUE	BLUE GREY	AQUA	GREEN	YELLOW	ORANGE	RED	VIOLET
HANDLING TINT											
5%											
10%											
20%											
30%											
40%											

**IMPORTANT NOTICE:** Colours shown here are representative of our tinted lenses and may not be 100% accurate. However, every effort has been made to ensure the closest match possible.

TINTS ARE AVAILABLE ON 3 DIFFERENT WATER CONTENTS: 58%, 67% AND 77%

58% = ● 67% = ● 77% = ●

TINT DENSITIES AVAILABLE	Handling Tint	5%	10%	20%	30%	40%
Colours						
Brown	●●●	●●●	●●●	●●●	●●	●
Green Brown	●●●	●●●	●●●	●●●	●●	●●
Grey	●●●	●●●	●●●	●●●	●●●	●●●
Blue Grey	●●●	●●●	●●●	●●●	●●●	●●
Blue		●●	●●	●●	●●	●●
Aqua	●●●	●●●	●●●	●●●	●●●	●●●
Green	●●●	●●●	●●●	●●●	●●●	●●●
Yellow				●●●		●●●
Orange				●●●		●●●
Red				●●●		●●●
Violet				●●●		●●●
Black	●●●					
Black Pupil	●●●					
Simulated Iris - Black	●●●					
Tinted Limbal Band*	●●					
Tinted Collarette*	●●					
Opaque	●●●					
Permanent Dot	●●●					

Warranty: We strongly recommend all tinted lenses are purchased with Warranty at the point of order.


\*Limbal bands required on 30% or 40% tint densities will only be supplied in black. Collarettes are not possible on tint densities of 30% and 40%.

It is not possible to have a Simulated Iris on a Toric design.


Please see reverse side for colour options.

TINTING OPTIONS


Tints are available on 3 different water contents: 58%, 67% and 77%. The material is a standard Hydrogel which is ideal for our range of tint colours and provides excellent comfort and stability.




**Type A**  
Iris area clear  
Pupil only tinted



**Type B**  
Whole lens tinted  
No separate pupil




**Type C**  
Iris area tinted  
Clear pupil




**Type D**  
Iris area tinted  
Pupil tinted separately

COLLARETTE AND LIMBAL BAND

UltraVision can supply a block colour Collarette and Limbal Band. The colours and densities available are as stated on the following page. Collarette and Limbal Bands are available on 67% and 77% material. Unless stated otherwise, the below rule will be applied to all limbal bands and collarettes:



Collarette



Limbal Band (1mm)

Handling Tint	= 5% Limbal Band/Collarette
5%	= 20% Limbal Band/Collarette
10%	= 30% Limbal Band/Collarette
20%	= 40% Limbal Band/Collarette
30%	= Black Limbal Band only
40%	= Black Limbal Band only


SIMULATED IRIDES

UltraVision supplies a range of tinted simulated irides for Prosthetic and Cosmetic contact lenses.


Detailed below are the different iris, collarette and limbal band print combinations available.

Please note Simulated Iris S1 will be supplied unless otherwise indicated.

All Simulated Irides are supplied in black, with opaque backing. All combinations can be supplied with a clear pupil.




Simulated Iris S1




Simulated Iris S1 with  
Limbal Band



Simulated Iris S2



Simulated Iris S2 with  
Collarette



Simulated Iris S2 with  
Collarette + Limbal  
Band

UltraVision CLPL Commerce Way, Leighton Buzzard, Bedfordshire LU74RW, United Kingdom

UK Orderline 0800 585115 (Freephone)

t +44 (0) 1525 381112 f +44 (0) 1525 370091

info@ultravision.co.uk ultravision.co.uk