

IGEL Pro

Simplify your life

A simple solution for complex patient needs - the **IGEL Pro** is a tailor-made soft contact lens ideal for high plus and minus prescriptions and delivering a tried and tested back surface toric design with cyl powers up to -11.00DC in 0.25DC steps.

The Empirical Fitting Tool ensures a simple fitting experience.

IGEL Pro is supplied in CONTAFLEX GM3 58, which is excellently tolerated, minimises depositing and ensures all day comfort.

Features

- Bi-curve, back surface toric-prism ballast
- Wide parameter range
- Empirical fitting tool available

Modality

Single lens
Daily wear
Up to 12 months.

Material

Contaflex GM3 58
58% water content

Warranty

First exchange - FOC
Second exchange - 50%
75% if lens unsuccessful

Base Curves (Rx Sphere)	6.40 to 10.00mm (0.10mm steps)
Base Curves (Rx Toric)	7.00 to 9.40mm (0.10mm steps)
Diameters	10.00 to 16.00mm (0.10mm steps)
Sphere Powers	-40.00 to +40.00 DS (0.25D steps)
Cyl Powers	-0.50 to -11.00DC (0.25D steps)
Axes	1° to 180° (1° steps)

Contact UltraVision to find out more...

ultravision
A MEMBER OF THE SEED GROUP

www.ultravision.co.uk

0800 585 115

info@ultravision.co.uk