



Which Solutions/Cleaners to use with UltraVision Lathe-cut SiH Contact Lenses

UltraVision recommends a rub and rinse system for cleaning any of our lathe-cut Silicone Hydrogel contact lenses, including:

- KeraSoft®3
- KeraSoft®IC (SiH)
- HydroWave®
- Aphakic SiH™
- Bandage SiH™
- High Myopic SiH™
- Paediatric SiH™

Below are guidelines on solutions deemed suitable for any of the above lenses:

Surfactant cleaner (Not alcohol based)

Disinfectant solutions, including:

- Focus Aqua
- OPTI-FREE® Express®
- OPTI-FREE® RepleniSH®
- ReNu Multiplus®
- Complete®
- Synergi
- MeniCare™ Soft

Peroxide Systems can be used if they are neutralised by either a platinum insert or tablet only:

- Oxysept
- Ever Clean®
- EasySEPT®

Enzyme Cleaners, including:

- Avizor Enzyme
- Amclair® (Abatron)

Please contact any one of our Customer Service team if you would like any further information.